



# Conjure the MAGIC of COMPRESSION PINNING

It's Not  
An Illusion



Imagine the lag screw effect  
with the ease of a k-wire  
- *in just one pin.*

Think of the time and effort  
you'll save.


The innovative thread design  
of the Orthofix Magic Pin  
auto-compresses bone  
fragments as the pin  
is advanced.

**We've made it so easy -  
it just seems like magic.**

# THE MAGIC PINS™ FRAGMENT FIXATION SYSTEM


Applications

## Small




- Middle and distal phalanges in the fingers and toes
- Fragments less than 3mm in diameter

## Medium



- Proximal phalanges in the fingers and toes
- Metacarpals and metatarsals
- Proximal and distal radius
- Small fragments in the humerus, tibia and fibula

## Large



- Humerus
- Olecranon
- Patella
- Tibia
- Fibula
- Metatarsals

Ordering Information

All small magic pins have:

- 120mm total length
- 1.5mm shaft diameter
- 1.2mm thread diameter

Small	Thread Length
W1207	7 mm
W1209	9 mm
W1211	11 mm
W1213	13 mm
W1215	15 mm
W1217	17 mm
W1219	19 mm
W1221	21 mm

All medium magic pins have:

- 120mm total length
- 2mm shaft diameter
- 1.6mm thread diameter

Medium	Thread Length
W1611	11 mm
W1613	13 mm
W1615	15 mm
W1617	17 mm
W1619	19 mm
W1621	21 mm
W1623	23 mm
W1625	25 mm
W1630	30 mm
W1635	35 mm
W1340	40 mm
W1645	45 mm
W1600	Washers*

All large magic pins have:

- 120mm total length
- 3mm shaft diameter
- 2.2mm thread diameter

Large	Thread Length
W2220	20 mm
W2225	25 mm
W2230	30 mm
W2235	35 mm
W2240	40 mm
W2245	45 mm
W2250	50 mm
W2255	55 mm
W2200	Washers*

*\*Note: Implants and washers are supplied in packs of 4*

W1001	Implant Extractor for Small and Medium implants	W1201	Sterilization box for small implants – Empty
W1002	Implant Extractor for Large implants	W1601	Sterilization box for medium implants – Empty
W1003	Implant Cutter	W2201	Sterilization box for large implants – Empty
W1000	A set of 4 empty Sterilization Boxes for small, medium and large implants and extractors and Allen wrench		
W1300	Complete set comprising: 40 small diameter Implants (5 Implants of each thread length); 60 medium diameter Implants (5 Implants of each thread length); 32 large diameter Implants (4 Implants of each thread length); 1 set of 4 Sterilization Boxes; 10 x W1600 Washers (medium); 10 x W2200 Washers (large); 1 x W1001 Implant Extractor (small/medium); 1 x W1002 Implant Extractor (large); 1 x 10.012 3mm Allen Wrench		



For more information, contact your local Orthofix representative, or call 1-800-BONEFIX

**ORTHOFIX®**  
www.orthofix.com