



Xtraordinary - packaged as a sterile kit in three time-saving configurations, ready for all long bone trauma applications.

Xtra light - as strong as conventional fixators, only lighter, for easier application and more patient comfort.

X-ray friendly - the first large fixator body to be radiolucent.

XCaliber - High Caliber Fixation from Orthofix.



XCaliber. **X**actly.

- Surprisingly Light Weight
- Excellent Stability
- Sterile Kit, Ready to Use



XCaliber Meta-Diaphyseal

XCaliber Meta-Diaphyseal

Stable Fixation of Metaphyseal and Diaphyseal Fractures

Multiple bone screw fixation possibilities with the same clamp.

Stable fixation (even in short bone segments) achieved through convergent and/or triangular screw positions.

Ball-joints allow correction up to 44° from straight position, in any direction, with fixator in place.

Ball-joints lock/unlock from either side. One instrument for fixator locking/unlocking.



Ideal in proximal and distal metaphysis.



Supplementary screw holder available for versatile screw placement.



Telescoping body for fracture reduction and compression/distraction. Dynamization for callus maturation.



Dynamic C/D Unit for compression /distraction and early cyclic movement for callus stimulation.

Case Study

37-year-old female, metaphyseal fracture (AO type C2).



Minimal internal fixation of the articular surface with cannulated screws. Metaphyseal bone graft with Norian SRS and application of the XCaliber Meta-Diaphyseal Fixator with three OsteoTite screws proximally and two distally.



Fixator removed at 140 days.



Case provided by: Dr. M. Manca, Hospital of Massa-Carrara, Italy.



Pre-assembled and Packaged Sterile
Ready to use, straight from the box.

Ordering Information

Part No.	Product Description
99-91215	XCaliber Meta-Diaphyseal Kit with compression/distraction unit, Sterile
99-91038	Supplementary Screw Holder Kit, Sterile
19200	Manipulation Forceps

XCaliber Fixator Is Available In Three Configurations.

Use with XCaliber Bone Screws.



Meta-Diaphyseal



Hybrid



Articulated Ankle

Your Distributor is:

www.orthofix.com

XC-0307(A)-PL-E0 05/04

ORTHOFIX
ALWAYS INNOVATING