

Prepared

with **PREFIX**



The Prefix™ fixator prepares you for any situation that requires rapid stabilization of complex fractures.

Prefix offers:

- Free pin placement in any desired plane
- Rapid, solid stabilization
- Simple reduction technique
- Radiolucency



Temporary joint bridging for ligamentotaxis



Rod Clip Clamp for quick connection of two rods




Screw Clip Clamp for free pin placement allows rapid assembly of desired frame construction

Multi-Screw Clamp for parallel screw positioning either in a T- or a straight clamp configuration



Radiolucent, lightweight, strong 14 mm diameter carbon fiber rod in four different lengths



Ordering Information

Screw Clip Clamp	92015
Rod Clip Clamp	92011
Transfixing Pin Clamp	92025
Transfixing Pin 4 mm, thread 5 x 50 mm	92050
Transfixing Pin 4 mm, thread 5 x 80 mm	92080
Rod 350 mm long	92350
Rod 275 mm long	92275
Rod 175 mm long	92175
Rod 125 mm long	92125
Allen Wrench 6 mm	10017
Multi-Screw Clamp	92030
Multi-Screw Clamp Guide	19930
PreFix Sterilization Box (empty)	92990

Indications

- Complex upper and lower extremity fractures
- Fractures in polytrauma
- Fractures in patients with pulmonary or cranial complications
- Open book pelvic fractures
- Temporary cross-leg fixation in plastic soft tissue reconstruction
- Temporary joint bridging

Your Distributor is:



www.orthofix.com