




MultiPlanar MiniRail™







MultiPlanar MiniRail™

The geared, M2 MultiPlanar MiniRail permits staged correction of complex joint deformities and non-surgical lengthening of musculotendinous contractures. With the M2 positioned over the center of rotation of the deformity, gradual correction in both the radio-ulnar and dorso-palmar planes and compression/distraction of either railed segment are possible.

Multiplanar Adjustment Micrometric Correction Mini Fixation




- Geared rails allow **gradual** positioning in both radio-ulnar or dorso-palmar planes
- Distraction and compression can be accomplished in either railed segment



Radial Club Hand Application

CODE	INDICATIONS
M511	<ul style="list-style-type: none"> ▪ Complex upper and lower extremity deformities ▪ Hand, wrist or foot contractures ▪ Indicated for paediatrics and adult patients

Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



www.orthofix.com

MR-0301-PL-E0 B 05/16

