


Excellent fracture visualisation, unobstructed view of radius and carpals



Radiolucent telescoping body

Smart and easy for optimal fracture x-ray visibility

Ball-joint design

Great stability and maneuverability to help anatomical reduction



19cm

13cm

Low profile design

Improved patient comfort



Compression/ Distraction unit

More accurate reduction

Sterile kit

Ready to use, with instrumentation, compression-distraction unit and standard bone screws



Ordering Information

Code No	Description
99-36501	RadioLucent Wrist Fixator Kit, sterile

For further ordering information see section "RadioLucent Wrist Fixator" in the Orthofix General Catalogue - See Quick Reference Guide PG 150

References:

1) Nelson D.L. : "The RadioLucent Wrist Fixator for Distal Radius Fractures" Orthofix External Fixation in Trauma and Orthopaedics, G. De Bastiani et al. Springer 2000, Ch. 17: 181-191.



Manufactured by:
 ORTHOFIX Srl
 Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
 Telephone +39 045 6719000, Fax +39 045 6719380



www.orthofix.com

PB RLW EO H 06/16

Distributed by:

