

Transverse Antegrade

Transverse Retrograde

Recon

Oblique

Short


Long

Antegrade

Retrograde


FEMORAL NAIL

Material Titanium
 Diameter 9-12 mm in 1 mm increments
 Length 275-475 mm in 25 mm increments
 Cannulation
 Ø 9 mm nail: no
 Ø 10-12 mm nail: yes
 Radius 2500 mm
 Locking bolts Ø 4.8 mm, unthreaded
 Recon/Oblique screws Ø 6.5 mm
 Condylar locking screws Ø 4.8 mm
 Color Red



TIBIAL NAIL

Material Titanium
 Diameter 8-11 mm in 1 mm increments
 Length 275-410 mm in 15 mm increments
 Cannulation
 Ø 8 mm nail: no
 Ø 9-11 mm nail: yes
 Bend 15°
 Locking bolts Ø 4.8 mm, unthreaded and threaded
 Ø 4.0 mm, unthreaded for Ø 8 mm nail
 Color Blue



SUPRACONDYLAR NAIL

Material Titanium
 Diameter 9-12 mm in 1 mm increments
 Length 130-375 mm in 25 mm increments
 Cannulation yes
 Bend 10°
 Radius 2500 mm
 Locking bolts Ø 4.8 mm, unthreaded
 Condylar locking screws Ø 4.8 mm
 Color Silver



HUMERAL NAIL

Proximal

Material Titanium
 Diameter 7-9 mm in 1 mm increments
 Length 150 mm
 Cannulation yes
 Bend 15°
 Locking bolts Ø 4.0 mm, unthreaded and threaded
 Color Gold

Standard

Material Titanium
 Diameter 7-9 mm in 1 mm increments
 Length 185-320 mm in 15 mm increments
 Cannulation yes
 Bend 15°
 Locking bolts Ø 4.0 mm, unthreaded and threaded
 Color Gold

